



When it's right,  
it just clicks

Amara

**PHILIPS**

RESPIRONICS

# Amara sleep apnea mask

Full face mask featuring an easy one "click" design

## Amara features just for you:

- Gel and silicone cushion options for your choice of comfort
- Lightweight and only four parts
- Easy to clean, easy to use
- Designed to fit up to 95% of users' faces
- All cushion sizes snap onto a single frame with a simple click



Fine-glide forehead adjuster with eight selectable positions

One "click" cushion assembly

Quick-fit headgear clips reduce the need for continual adjustment

Lightweight and only four parts

360° rotating elbow with quick release swivel



Learn more about Amara at [www.philips.us/youramara](http://www.philips.us/youramara)

**Philips Healthcare is part of  
Royal Philips**

How to reach us

[www.philips.com/healthcare](http://www.philips.com/healthcare)

Philips Respironics  
1010 Murry Ridge Lane  
Murrysville, PA 15668

Respironics and Amara are trademarks of Respironics, Inc. and its affiliates.

Please visit [www.philips.us/youramara](http://www.philips.us/youramara)



© 2015 Koninklijke Philips N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from use of this publication.  
CAUTION: U.S. federal law restricts these devices to sale by or on the order of a physician.

RRDPGH MG 12/10/14 MCI 4106508 PN 1110439